




BruTab_{6S}[®]

**Effervescent Disinfectant
& Sanitizer Tablet**



How NaDCC Works

How it differs from traditional Hypochlorite Bleach

The active ingredient in **BruTab 6S** is Sodium dichloro-s-triazinetriene (NaDCC), which produces an Available Chlorine solution that is very effective as a disinfecting and sanitizing agent against a broad spectrum of micro-organisms.

While **BruTab 6S** does provide chlorine, it is not a hypochlorite as is traditional bleach. There are significant differences that need to be understood to differentiate the two product types.

Sodium hypochlorite in bleach is stabilized with caustics and as a result, it has a pH of 11 or higher. NaDCC on the other hand has a neutral pH of ~6.5 when dissolved in water making NaDCC less harmful to surfaces and skin than acidic and alkaline disinfectants.

Once in solution, NaDCC releases approximately 50% of its total chlorine content as Free Available Chlorine (FAC) which is the active disinfection agent. As the FAC is consumed during the disinfection process, the NaDCC continues to release chlorine maintaining the 50/50 equilibrium in solution for longer lasting disinfection power than bleach. Sodium hypochlorite releases all its chlorine content immediately and once consumed there is no replenishment making it less effective. This means that **BruTab 6S** is not inactivated by dirt/cloths/organic matter as easily as bleach.

FAC exists in two forms, Hypochlorous acid (HOCl) found in **BruTab 6S** solutions and Hypochlorite ion (OCl-) found in a bleach solution. Studies have shown that Hypochlorous acid has 4X (four times) more disinfection power than the hypochlorite ion. HOCl is very similar to the water molecule allowing it to easily penetrate through the negatively charged cell wall. Once the HOCl enters the micro-organism, it destroys the nucleus of the cell completing the disinfection process. The Hypochlorous acid found in bleach is used up very quickly, it doesn't penetrate the cell as easily and solutions must be replaced more often, especially in the presence of organic loads.

These differences lead to the unique advantages of **BruTab 6S**:

- Delivers more potent, longer lasting disinfection power in the form of Hypochlorous acid – **Strong and Cost Effective**
- Has long lasting, available chlorine in reserve – **Stable** solution that can be stored for 3 days in a closed container. Bleach becomes inactive after a day.
- Its tablet form is **Stable** for 3 years from date of manufacture. Bleach concentrates have a 6 month shelf life – **Stable/Sustainable** product
- Neutral pH is less damaging to surfaces and equipment than acidic and alkaline disinfectants.



BruTab_{6S}[®]

1

STRONG

One product for use throughout an entire facility

US EPA Registered Cleaner-Disinfectant with efficacy against *C. difficile* spores, TB, and other microorganisms in 4 minutes. Kills Norovirus, Hepatitis A Virus, Hepatitis B Virus, Hepatitis C Virus and HIV-1 in 1 minute. Hospital Disinfectant. Complies with surface disinfection requirements of OSHA Bloodborne Pathogens Standard. Sanitizer Claim for Food Service Applications; Kills *Canine Parvovirus*.

Dilution and contact times vary by organism.

2

SURFACE/USER FRIENDLY

Neutral pH is less damaging to surfaces and equipment than acidic and alkaline disinfectants.

Neutral pH similar to skin. Will not dull floor finishes. Safe on most colorfast fabrics.


** Always test small area of clothing for color fastness before using.*

3

STABLE

Economical — lasts longer resulting in less waste

Longer shelf life in solution than bleach – 3 days compared to 1 day for bleach. Longer shelf life in tablet form than bleach – years compared to months for bleach. Continues working in the presence of organic load (i.e. blood and dirt). Sodium dichloro-s-triazinetriene retains killing power due to a 50 / 50 chemical equilibrium that continues to generate hypochlorous acid to replace that which is being used up in the process of destroying micro-organisms or contact with organic loads. This means that it is not inactivated by dirt/cloths/organic matter as easily as chlorine is depleted from bleach.



4

SUSTAINABLE

Cost-savings in multiple ways: Storage, Shipping, Handling, Waste Minimization

Small tablet size: Reduces SKUs – less warehouse space required to stock product compared to bleach and other liquid disinfectants. Reduces shipping costs. Non-hazardous shipping. Reduces packaging waste.

5

SIMPLE

Ease of training and usage

Exact dosage every time delivers accurate strength solution. Eliminates “measure and pour” guesswork. No dispensing equipment required.

6

SMELLS CLEAN

Leaves areas smelling clean and disinfected

Ideal for all areas in the facility including patient rooms, restrooms and shower rooms.

13.1g Tablets



PACKAGING INFORMATION

Availability

- Product No. 161021
 • 256 Tablet Tub/2 (8N)

BruTab 6S® is a US EPA registered broad spectrum disinfectant, virucide and food contact surface sanitizer as has been demonstrated by its performance in tests that are prescribed and regulated by the federal government under the Federal Insecticide, Fungicide, and Rodenticide Act (FIFRA).

DILUTION CHART

| Tablet Size | 13.1 g | |
|--|---------|------------------|
| Solution ppm (mg/L) Available Chlorine | Tablets | Gallons of Water |
| 100 | 1 | 10 |
| 538 | 1 | 2 |
| 1076 | 1 | 1 |
| 2153 | 2 | 1 |
| 4306 | 4 | 1 |
| 5382 | 5 | 1 |

Use-dilution labels available for prepared solution.

EPA Reg. No. 71847-6-106
 EPA 27 DEC 18

Refer to label for directions for use, claims and other organisms.

TESTING SUMMARY:

| Human Microorganisms for Disinfection, Germicidal, Virucidal, and Fungicidal Claims | | ATCC and/or Strain Number and/or Surrogate organism | Disease/Effect | 1 minute minimum contact time with organic soil load # tablets (ppm Solution) | 4 minute minimum contact time with organic soil load # tablets (ppm Solution) | 10 minute minimum contact time # tablets (ppm Solution) |
|---|----------|---|------------------------|---|---|---|
| <i>Clostridium difficile</i> spores | Spores | ATCC 43598 | Colitis | | 4 per gal. (4306 ppm) | 2 per gal. (2153 ppm) |
| <i>Mycobacterium bovis</i> (TB) | Bacteria | ATCC 35743 | Tuberculosis (TB) | | 5 per gal. (5382 ppm) | |
| <i>Acinetobacter baumannii</i> | Bacteria | ATCC BAA-1709 | Wound infections etc. | | 4 per gal. (4306 ppm) | |
| Vancomycin Resistant <i>Enterococcus faecalis</i> (VRE) | Bacteria | ATCC 51575 | Enteritis etc. | | 4 per gal. (4306 ppm) | 1 per gal. (1076 ppm) |
| <i>Escherichia coli</i> 0157:H7 | Bacteria | ATCC 35150 | Food poisoning | | | 1 per gal. (1076 ppm) |
| <i>Klebsiella pneumoniae</i> | Bacteria | ATCC 4352 | Pneumonia | | | 1 per gal. (1076 ppm) |
| Carbapenem resistant <i>Klebsiella pneumoniae</i> | Bacteria | ATCC BAA-1705 | Pneumonia | | 4 per gal. (4306 ppm) | |
| <i>Pseudomonas aeruginosa</i> | Bacteria | ATCC 15442 | Septicemia | | 4 per gal. (4306 ppm) | 1 per gal. (1076 ppm) |
| <i>Salmonella enterica</i> | Bacteria | ATCC 10708 | Food poisoning | | 4 per gal. (4306 ppm) | 1 per gal. (1076 ppm) |
| <i>Staphylococcus aureus</i> | Bacteria | ATCC 6538 | Wound infections etc. | | 4 per gal. (4306 ppm) | 1 per gal. (1076 ppm) |
| <i>Staphylococcus aureus</i> MRSA & GRSA | Bacteria | ATCC 33592 | Wound infections etc. | | 4 per gal. (4306 ppm) | 1 per gal. (1076 ppm) |
| <i>Staphylococcus epidermis</i> | Bacteria | ATCC 51624 | Wound infections, etc. | | | 1 per gal. (1076 ppm) |
| <i>Streptococcus pneumoniae</i> | Bacteria | ATCC 6305 | Pneumonia | | 4 per gal. (4306 ppm) | |

TESTING SUMMARY:

| Human Microorganisms for Disinfection, Germicidal, Virucidal, and Fungicidal Claims | | ATCC and/or Strain Number and/or Surrogate organism | Disease/Effect | 1 minute minimum contact time with organic soil load # tablets (ppm Solution) | 4 minute minimum contact time with organic soil load # tablets (ppm Solution) | 10 minute minimum contact time # tablets (ppm Solution) |
|---|-----------------------|---|------------------------------|---|---|---|
| Coxsackievirus | Virus (Non-Enveloped) | ATCC VR-30 | Hand, foot and mouth disease | 4 per gal. (4306 ppm) | | |
| Hepatitis A Virus | Virus (Non-Enveloped) | Strain HM-175/18f | Hepatitis A | 4 per gal. (4306 ppm) | | 1 per gal. (1076 ppm) |
| Norovirus | Virus (Non-Enveloped) | ATCC VR-782 | Gastroenteritis | 2 per gal. (2153 ppm) | | |
| Poliovirus Type 1 | Virus (Non-Enveloped) | ATCC VR-1000 | Polio | | | 1 per gal. (1076 ppm) |
| Rhinovirus (Type 14) | Virus (Non-Enveloped) | ATCC VR284 | Common cold | | | 1 per gal. (1076 ppm) |
| Avian Influenza A Virus (H5N1) | Virus (Enveloped) | CDC #2006719965 | Flu | 4 per gal. (4306 ppm) | | |
| Hepatitis B Virus | Virus (Enveloped) | Duck Hepatitis B (DHBV) | Hepatitis B | 4 per gal. (4306 ppm) | | 1 per gal. (1076 ppm) |
| Hepatitis C Virus | Virus (Enveloped) | Bovine Viral Diarrhea Virus Strain NADL | Hepatitis C | 4 per gal. (4306 ppm) | | |
| Herpes Simplex Virus Type 1 | Virus (Enveloped) | ATCC VR-733 | Herpes | | | 1 per gal. (1076 ppm) |
| Human Immunodeficiency Virus Type 1 | Virus (Enveloped) | Strain IIIB | AIDS | 4 per gal. (4306 ppm) | | 1 per gal. (1076 ppm) |
| Influenza Virus (H1N1) | Virus (Enveloped) | ATCC VR-99 | Swine flu | | | 1 per gal. (1076 ppm) |
| Respiratory syncytial virus | Virus (Enveloped) | ATCC VR-26 | Common cold | | | 1 per 2 gal. (538 ppm) |
| Aspergillus fumigatus | Fungi | ATCC 36607 | Respiratory infections | 4 per gal. (4306 ppm) | | |
| Candida albicans | Fungi | ATCC 10231 | Thrush & Yeast Infections | | 4 per gal. (4306 ppm) | |
| Trichophyton interdigitale | Fungi | ATCC 9533 | Athlete's foot | | | 1 per gal. (1076 ppm) |

| Animal Microorganisms | | ATCC and/or Strain Number | Disease/Effect | 1 minute contact time with organic soil load # tablets (ppm Solution) | 4 minute contact time with organic soil load # tablets (ppm Solution) | 10 minute contact time # tablets (ppm Solution) |
|-------------------------|-----------------------|---------------------------|--------------------|---|---|---|
| Canine Parvovirus | Virus (Non-Enveloped) | ATCC VR-2017 | Parvovirus disease | | | 1 per gal. (1076 ppm) |
| Feline Calicivirus | Virus (Non-Enveloped) | ATCC VR-782 | Gastroenteritis | 2 per gal. (2153 ppm) | | 1 per gal. (1076 ppm) |
| Canine Distemper Virus | Virus (Enveloped) | ATCC VR-128 | Canine distemper | | | 1 per gal. (1076 ppm) |
| Newcastle Disease Virus | Virus (Enveloped) | ATCC VR-180 | Newcastle disease | | | 1 per gal. (1076 ppm) |
| Pseudorabies Virus | Virus (Enveloped) | ATCC VR-135 | Aujesky's disease | | | 1 per gal. (1076 ppm) |

See page 10 for complete list

| Food Contact Surface Sanitizer When used at 100 ppm solution, applied as outlined under Sanitizer Directions, BruTab 6S is an effective food contact surface sanitizer | | ATCC and/or Strain Number | Disease/Effect | 1 minute contact time # tablets (ppm Solution) |
|---|----------|---------------------------|-----------------------|--|
| Salmonella enterica | Bacteria | ATCC 6539 | Food poisoning | 1 per 10 gal. (100 ppm) |
| Staphylococcus aureus | Bacteria | ATCC 6538 | Wound infections etc. | 1 per 10 gal. (100 ppm) |

3.3g Tablets



PACKAGING INFORMATION

Availability

Product No. 161021
 • 200 Tablet Tub/6 (8G)

BruTab 6S® is a US EPA registered broad spectrum disinfectant, virucide and food contact surface sanitizer as has been demonstrated by its performance in tests that are prescribed and regulated by the federal government under the Federal Insecticide, Fungicide, and Rodenticide Act (FIFRA).

DILUTION CHART

| Tablet Size | 3.3g | |
|--|---------|-----------------|
| Solution ppm (mg/L) Available Chlorine | Tablets | Quarts of Water |
| 100 | 1 | 10 |
| 538 | 1 | 2 |
| 1076 | 1 | 1 |
| 2153 | 2 | 1 |
| 4306 | 4 | 1 |
| 5382 | 5 | 1 |

Use-dilution labels available for prepared solution.

EPA Reg. No. 71847-6-106
 EPA 27 DEC 18

Refer to label for directions for use, claims and other organisms.

TESTING SUMMARY:

| Human Microorganisms for Disinfection, Germicidal, Virucidal, and Fungicidal Claims | | ATCC and/or Strain Number and/or Surrogate organism | Disease/Effect | 1 minute minimum contact time with organic soil load # tablets (ppm Solution) | 4 minute minimum contact time with organic soil load # tablets (ppm Solution) | 10 minute minimum contact time # tablets (ppm Solution) |
|---|----------|---|------------------------|---|---|---|
| <i>Clostridium difficile</i> spores | Spores | ATCC 43598 | Colitis | | 4 per quart (4306 ppm) | 2 per quart (2153 ppm) |
| <i>Mycobacterium bovis</i> (TB) | Bacteria | ATCC 35743 | Tuberculosis (TB) | | 5 per quart (5382 ppm) | |
| <i>Acinetobacter baumannii</i> | Bacteria | ATCC BAA-1709 | Wound infections etc. | | 4 per quart (4306 ppm) | |
| Vancomycin Resistant <i>Enterococcus faecalis</i> (VRE) | Bacteria | ATCC 51575 | Enteritis etc. | | 4 per quart (4306 ppm) | 1 per quart (1076 ppm) |
| <i>Escherichia coli</i> 0157:H7 | Bacteria | ATCC 35150 | Food poisoning | | | 1 per quart (1076 ppm) |
| <i>Klebsiella pneumoniae</i> | Bacteria | ATCC 4352 | Pneumonia | | | 1 per quart (1076 ppm) |
| Carbapenem resistant <i>Klebsiella pneumoniae</i> | Bacteria | ATCC BAA-1705 | Pneumonia | | 4 per quart (4306 ppm) | |
| <i>Pseudomonas aeruginosa</i> | Bacteria | ATCC 15442 | Septicemia | | 4 per quart (4306 ppm) | 1 per quart (1076 ppm) |
| <i>Salmonella enterica</i> | Bacteria | ATCC 10708 | Food poisoning | | 4 per quart (4306 ppm) | 1 per quart (1076 ppm) |
| <i>Staphylococcus aureus</i> | Bacteria | ATCC 6538 | Wound infections etc. | | 4 per quart (4306 ppm) | 1 per quart (1076 ppm) |
| <i>Staphylococcus aureus</i> MRSA & GRSA | Bacteria | ATCC 33592 | Wound infections etc. | | 4 per quart (4306 ppm) | 1 per quart (1076 ppm) |
| <i>Staphylococcus epidermis</i> | Bacteria | ATCC 51624 | Wound infections, etc. | | | 1 per quart (1076 ppm) |
| <i>Streptococcus pneumoniae</i> | Bacteria | ATCC 6305 | Pneumonia | | 4 per quart (4306 ppm) | |

TESTING SUMMARY:

| Human Microorganisms for Disinfection, Germicidal, Virucidal, and Fungicidal Claims | | ATCC and/or Strain Number and/or Surrogate organism | Disease/Effect | 1 minute minimum contact time with organic soil load # tablets (ppm Solution) | 4 minute minimum contact time with organic soil load # tablets (ppm Solution) | 10 minute minimum contact time # tablets (ppm Solution) |
|---|-----------------------|---|------------------------------|---|---|---|
| Coxsackievirus | Virus (Non-Enveloped) | ATCC VR-30 | Hand, foot and mouth disease | 4 per quart (4306 ppm) | | |
| Hepatitis A Virus | Virus (Non-Enveloped) | Strain HM-175/18f | Hepatitis A | 4 per quart (4306 ppm) | | 1 per quart (1076 ppm) |
| Norovirus | Virus (Non-Enveloped) | ATCC VR-782 | Gastroenteritis | 2 per quart (2153 ppm) | | |
| Poliovirus Type 1 | Virus (Non-Enveloped) | ATCC VR-1000 | Polio | | | 1 per quart (1076 ppm) |
| Rhinovirus (Type 14) | Virus (Non-Enveloped) | ATCC VR284 | Common cold | | | 1 per quart (1076 ppm) |
| Avian Influenza A Virus (H5N1) | Virus (Enveloped) | CDC #2006719965 | Flu | 4 per quart (4306 ppm) | | |
| Hepatitis B Virus | Virus (Enveloped) | Duck Hepatitis B (DHBV) | Hepatitis B | 4 per quart (4306 ppm) | | 1 per quart (1076 ppm) |
| Hepatitis C Virus | Virus (Enveloped) | Bovine Viral Diarrhea Virus Strain NADL | Hepatitis C | 4 per quart (4306 ppm) | | |
| Herpes Simplex Virus Type 1 | Virus (Enveloped) | ATCC VR-733 | Herpes | | | 1 per quart (1076 ppm) |
| Human Immunodeficiency Virus Type 1 | Virus (Enveloped) | Strain IIIB | AIDS | 4 per quart (4306 ppm) | | 1 per quart (1076 ppm) |
| Influenza Virus (H1N1) | Virus (Enveloped) | ATCC VR-99 | Swine flu | | | 1 per quart (1076ppm) |
| Respiratory syncytial virus | Virus (Enveloped) | ATCC VR-26 | Common cold | | | 1 per 2 quarts (538 ppm) |
| Aspergillus fumigatus | Fungi | ATCC 36607 | Respiratory infections | 4 per quart (4306 ppm) | | |
| Candida albicans | Fungi | ATCC 10231 | Thrush & Yeast Infections | | 4 per quart (4306 ppm) | |
| Trichophyton interdigitale | Fungi | ATCC 9533 | Athlete's foot | | | 1 per quart (1076 ppm) |

| Animal Microorganisms | | ATCC and/or Strain Number | Disease/Effect | 1 minute contact time with organic soil load # tablets (ppm Solution) | 4 minute contact time with organic soil load # tablets (ppm Solution) | 10 minute contact time # tablets (ppm Solution) |
|-------------------------|-----------------------|---------------------------|--------------------|---|---|---|
| Canine Parvovirus | Virus (Non-Enveloped) | ATCC VR-2017 | Parvovirus disease | | | 1 per quart (1076 ppm) |
| Feline Calicivirus | Virus (Non-Enveloped) | ATCC VR-782 | Gastroenteritis | 2 per quart (2153 ppm) | | 1 per quart (1076 ppm) |
| Canine Distemper Virus | Virus (Enveloped) | ATCC VR-128 | Canine distemper | | | 1 per quart (1076 ppm) |
| Newcastle Disease Virus | Virus (Enveloped) | ATCC VR-180 | Newcastle disease | | | 1 per quart (1076 ppm) |
| Pseudorabies Virus | Virus (Enveloped) | ATCC VR-135 | Aujesky's disease | | | 1 per quart (1076 ppm) |

See page 10 for complete list

| Food Contact Surface Sanitizer When used at 100 ppm solution, applied as outlined under Sanitizer Directions, BruTab 6S is an effective food contact surface sanitizer | | ATCC and/or Strain Number | Disease/Effect | 1 minute contact time # tablets (ppm Solution) |
|---|----------|---------------------------|-----------------------|--|
| Salmonella enterica | Bacteria | ATCC 6539 | Food poisoning | 1 per 10 quart (100 ppm) |
| Staphylococcus aureus | Bacteria | ATCC 6538 | Wound infections etc. | 1 per 10 quart (100 ppm) |



BruMop

Bucketless Handle Mop System

Easy to use bucketless handle mop system with 32 fl oz reservoir to use for mopping floors with BruTab 6S. Use one 13.1g tablet in 32 fl oz reservoir to prepare 4306 ppm solution. Use one 13.1g tablet with 24 fl oz water to prepare 5382 ppm solution for 4 minute *Mycobacterium bovis* (TB) claim.


Availability

Product No. ATWK18OR

- One orange applicator with 32 fl oz reservoir. (does not ship with a microfiber pad)

Product No. ATWR18OR-CS

- 18" orange microfiber flat mop (CS/12)



BruTab 6S Transport Pack

Ten empty containers along with a sheet of customer affixed labels for distributing tablets to different areas within a facility.

Availability

Product No. PDU102



ANIMAL PATHOGENS

When used at ppm solution listed below, applied as outlined under Animal Premises Disinfection/ Virucidal Directions, BruTab 6S is effective against the following animal pathogens with the corresponding contact time:

10 minute contact time (1076 ppm Solution)

Canine Parvovirus
 Newcastle Disease Virus
 Pseudorabies Virus
 Canine Distemper Virus
 Feline Calicivirus
 Actinobacillus pleuropneumoniae
 Bordetella bronchiseptica (Rhinitis)
 Brachyspira (Treponema/Serpulina)
 Gumboro disease Virus
 Herpes Simplex Virus Type 1
 Hyodysenteriae (Swine Dysentery)
 Infectious Canine hepatitis Virus
 Porcine epidemic diarrhea Virus
 Porcine parvovirus
 Runting & Stunting Virus (tenosynovitis)
 Streptococcus uberis
 Teschen/Talfan disease Virus
 Respiratory syncytial virus (538 ppm)

1 minute contact time (ppm Solution)

Feline Calicivirus (2153 ppm)
 Avian influenza virus (H5N1) (4306 ppm)
 Bovine Viral Diarrhea Virus (4306 ppm)

NOTE: Only approved for use against Canine Parvovirus, Newcastle Disease Virus, Pseudorabies, Canine Distemper Virus and Feline Calicivirus in the State of California.

PHYSICAL & CHEMICAL SPECIFICATIONS

| | |
|---|-----------------|
| Active Ingredient: Sodium dichloro-s-triazinetriene | 48.21%* |
| Working pH | 6.0 +/- 0.5 |
| Color | Clear |
| Odor | Slight Chlorine |
| OSHA GHS Rating In-Use | Non-Hazardous |

** Equivalent to 31.10% active chlorine by tablet weight.
Refer to dilution chart for Available Chlorine concentrations*

STABILITY DATA

A stability study was conducted and found that solutions made up of strengths from 100 - 5382 ppm active chlorine, retained the required chlorine activity after storage for 3 days in a closed container at room temperature out of direct sunlight.

General Solution Preparation: Prepare a fresh solution twice weekly (every 3 days) when using closed containers (spray bottles). Prepare a fresh solution daily when using open containers (buckets) or if solution becomes diluted.



MATERIAL SUBSTRATE COMPATIBILITY

Sodium dichloro-s-triazinetriene tablets dissolved in water produce a solution of hypochlorous acid.

The following chart shows the compatibility of a variety of materials with solutions up to 5,000 mg/L of available chlorine.

| Metals | Compatibility |
|--------------|---------------|
| SS 304 | A |
| SS 316 | A |
| Aluminum | B |
| Brass | B |
| Bronze | B |
| Carbon Steel | C |
| Cast iron | C |
| Hasteloy C® | A |
| Titanium | A |

The following chart shows the compatibility of a variety of materials with solutions up to 2,000 mg/L of available chlorine.

| Plastics | Compatibility | Elastomers | Compatibility |
|---------------|---------------|------------------|---------------|
| ABS | A | Nitrile (Buna N) | A |
| CPVC | A | EPDM | A |
| Hytrel® | A | Hypalon® | A |
| HDPE | A | Kel-F® | A |
| LDPE | A | Santoprene® | A |
| Noryl® | A | Silicone | B |
| Polycarbonate | A | Tygon® | A |
| Polypropylene | A | Viton® | A |
| PPS | A | Nonmetals | Compatibility |
| PTFE | A | Carbon graphite | A |
| PVC | A | Ceramic A1203 | A |
| PVDF | A | Ceramic magnet | A |

Explanation of Ratings — Chemical Effect

A = Excellent.

B = Good – Minor Effect, slight corrosion or discoloration.

C = Fair – Moderate Effect, OK for short term use. Not recommended for continuous use. Some pitting may occur.

D = Severe Effect, not recommended for any use.





317.923.3211
WWW.BHCINC.COM
ISO 9001:2015 CERTIFIED

*BruLin and BruTab 6S are registered
trademarks of BHC, Inc.*

Distributed by: